Levita

compressor nebulizer

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
INSTRUCTIONS FOR USE MANUAL

Aerosol Therapy Device

Congratulations on your purchase and thank you for your confidence in us. Our goal is to fully satisfy our customers by offering them cutting-edge products for treating ailments of the respiratory system. **Please read these instructions carefully and keep them handy for future reference. Use this device only as described in this manual. This is a home medical device for the nebulization and administration of medication prescribed or suggested by your doctor.**

THIS DEVICE COMES EQUIPPED WITH:


A – Aerosol therapy unit (main unit)

- A1 - Switch
- A2 - Air outlet port
- A3 - Air filter
- A4 - Nebulizer port
- A5 - Carrying handle
- A6 - Power cord compartment
- A7 - Socket for power cord
- A8 - Power cord
- A9 - Cover
- A10 - Accessory compartment

B - Connection tube (main unit / nebulizer)

C - Accessories

- C1 – RF8 Plus nebulizer
 - C1.1 – Lower part
 - C1.2 – Upper part with valve and speed selector
 - C1.2.1 – Upper part
 - C1.2.2 – Inspiratory valve system
 - C1.2.2 a – Inspiratory valve
 - C1.2.2 b – Valve retainer
 - C1.2.3 – Speed selector
- C2 – Nasal prong
- C3 – Mouthpiece with valve
 - C3.1 Expiratory valve
- C4 – Adult mask
- C5 – Child mask
- C6 – Nebulization manual control

IMPORTANT SAFEGUARDS

- Before using the device for the first time and periodically during the product's life, check the integrity of its structure and of the power cable, making sure that there is no damage; if it is damaged, do not insert the plug and take the product immediately to an authorised technical assistance center or to your local dealer.
- The average lifespan of accessories is 1 year; however, we recommend you replace the nebulizer every 6 months in intensive use (or earlier, if the nebulizer is blocked) to always ensure maximum therapeutic effect.
- Children and disabled persons should always use the unit under strict supervision of an adult who has read this manual.
- Several parts of the device are small enough to be swallowed by children; therefore keep the device out of reach of children.
- Do not use the tubes and the cords of the device for purposes other than those specified as they may cause strangulation hazard; utmost attention must be paid to children and persons with special needs; usually, such persons are not capable of assessing the hazard correctly.
- The unit should not be used in the presence of anesthetic mixture inflammable with air, oxygen or nitrous oxide.
- Keep the cord away from heated surfaces.
- Keep the power supply cord away from animals (for example rodents) as they may damage the insulation of the power supply cord.
- Do not handle the unit with wet hands. Do not use the device in humid environments (e.g. while taking a bath or shower). Do not immerse the device in water; should this accidentally occur, unplug the device immediately. Do not remove or touch the device immersed in water before pulling out the plug. Take the device immediately to an authorized service center or to your local dealer.
- Only use the device in dust-free environments; otherwise, therapy may be compromised.
- The casing of this unit is not protected against the penetration of liquids. Do not wash the device under running water or submerge in water. Keep away from splashing water or other liquids.
- Do not expose the device to extreme temperatures.
- Do not put the device near heat sources, direct sunlight, or in hot environments.
- Do not obstruct or insert any objects in the filter or in its housing in the device.
- Do not obstruct the air vents placed on both sides of the unit.
- During use, always place the unit on a firm surface clear of obstacles.
- Check that there is no material obstructing the air vents before each use.
- Do not place any objects inside the air vents.
- Repairs must be performed by authorized personnel only. Any unauthorized repair will void the warranty and may pose a safety hazard for the user.

WARNING: Do not modify this equipment without authorization of the manufacturer

- The average operational life of the compressor is 700 hours.
- The Manufacturer, the Retailer and the Importer are responsible for the safety, reliability and performance of the device only if: a) the device is used in accordance with the operating instructions b) the electrical system of the premises where the device is used is in compliance with current laws.
- Interactions: The materials used that come into contact with medications have been tested with a wide range of medications. However, due to the variety and continuous evolution of medication, the possibility of a chemical interaction cannot be excluded.

We recommend that once the medication has been opened you should use it as soon as possible and avoid prolonged contact of the medication with the nebulizer.

- The manufacturer should be contacted for reporting problems and/or unexpected events related to device operation.

INSTRUCTIONS FOR USE

Before each use, the nebulizer and accessories should be cleaned according to the instructions in the section "CLEANING, SANITIZATION AND DISINFECTION". A personal use of the nebulizer and its accessories is recommended in order to avoid any risks of infection.

During application, it is advisable to suitably protect yourself from any dripping.

1. Plug the power cord (A8) into the unit socket (A7) and plug the power cord into a power socket corresponding to the voltage of the unit. The power cord must be positioned so that it is not difficult to make the disconnection from the mains.

2. Assemble the nebulizer upper part, as shown in the assembly diagram section C1.2.

Warning: The arrow shown on the valve retainer (C1.2.2) should be aligned with the accessory connection port (see picture). This will ensure an accurate fit into the upper part (C1.2.1). Pour the medication prescribed by the doctor into the lower part (C1.1). Insert the upper part (C1.2) into the lower part (C1.1), close the nebulizer by turning the upper part (C1.2) clockwise.


3. Connect accessories as indicated in the "Assembly diagram" on the cover.

Sit comfortably holding the nebulizer in your hand, place the mouthpiece into your mouth, alternatively you can use the nose piece or mask. Should you opt for the mask accessory, place it on your face as shown in the picture (with or without using the elastic).

4. Switch on the unit (A1) Gently take a deep breath; we recommend holding your breath for an instant to allow the inhaled aerosol droplets to deposit. Then exhale slowly.

5. When the treatment is finished, switch off and unplug the unit.

WARNING. If after the therapy session moisture accumulates in the tube (B), remove the tube from the nebulizer and dry it using the compressor air flow; this operation will prevent mold from forming inside the tube.

OPERATION MODES OF THE RF8 PLUS NEBULIZER WITH SPEED SELECTOR AND VALVE SYSTEM

This nebulizer chamber is professional, fast and suitable to administer all medication.


To make the inhalation therapy faster, turn the speed selector (C1.2.3) clockwise.

For more effective inhalation therapy, turn the speed selector (C1.2.3) anticlockwise.

The valve system in the nebulizer and in the mouthpiece will allow optimization of the medication intake, limiting the medication dispersion into the environment


SoftTouch Maske

SoftTouch masks have an outer edge made of soft biocompatible material that ensures excellent adherence to the face, and is also equipped with an **innovative Dispersion Limiting**. These innovations allow increased sedimentation of medication in the lungs and reduced dispersion into the air.


During inhaling
the vent bends inwards.

During exhaling
the vent bends outwards.

USE OF NEBULIZATION MANUAL CONTROL

The nebulisation manual control (C6) can be used to reduce dispersion of the medication into the surrounding environment and reduces wastage of the medication. It should not be used if continuous nebulisation of the medication is required, which may be more suited for children and people with reduced capabilities.

To start nebulizing close with a finger the hole of the nebulization manual control (C6) and breathe in gently; we recommend to hold your breath for a moment so that the inhaled aerosol droplets can be deposited,

meanwhile, to disable nebulizing, remove your finger from the hole of the nebulization manual control to avoid the waste of drug, optimizing its acquisition. Then exhale slowly.

CLEANING, SANITIZATION AND DISINFECTION

Turn off the unit before any cleaning and disconnect the power cord from the mains.

DEVICE AND OUTER SURFACE OF THE TUBE

Only use a cloth dampened with antibacterial detergent (non-abrasive and free of solvents of any kind).

ACCESSORIES

We recommend a personal use of both the nebulizer and accessories to avoid any risk of infection or contagion. Open the nebulizer by turning the upper part (C1.2) anticlockwise, then dismantle it as shown in the "Connection diagram" sect. C1. Then proceed according to the following instructions.

SANITIZATION

Before and after each use, sanitize the nebulizer and accessories according to one of the methods described below.

method A: Sanitize accessories **C1.1-C1.2.1-C1.2.2-C1.2.3- C2-C3-C4-C5-C6** under potable hot water (approximately 40°C) with a gentle, non abrasive dish detergent.

method B: Sanitize accessories **C1.1-C1.2.1-C1.2.2-C1.2.3- C2-C3-C4-C5-C6** in the dishwasher with a hot cycle.

method C: Sanitize accessories **C1.1-C1.2.1-C1.2.2-C1.2.3- C2-C3-C4-C5-C6** by immersing them in a solution of 50% water and 50% white vinegar, then rinse thoroughly under potable hot water (approximately 40°C).

If you want to also perform DISINFECTION, proceed to the DISINFECTION paragraph.

After sanitization shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).

DISINFECTION

After sanitizing the nebulizer and the accessories, disinfect them choosing one of the methods described below.

method A: Accessories **C1.1-C1.2.1-C1.2.2-C1.2.3- C2-C3-C4-C5-C6** can be disinfected.

The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite), specific for disinfecting, which is available in any pharmacy.

Procedure:

- Fill a container, suitable to contain all the individual components to be disinfected, with a solution of drinking water and disinfectant, while respecting the proportions indicated on the packaging of the disinfectant itself.

- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the packaging of the disinfectant, and associated with the concentration chosen to prepare the solution.

- Remove the components now disinfected and rinse thoroughly with lukewarm drinking water.

- Dispose of the solution as recommended by the disinfectant manufacturer.


method B: Disinfect the accessories **C1.1-C1.2.1-C1.2.2-C1.2.3-C2-C5-C6** by boiling them in water for 10 minutes; use demineralised or distilled water to prevent calcium deposits.

method C: Disinfect the accessories **C1.1-C1.2.1-C1.2.2-C1.2.3-C2-C5-C6** with a hot steam sterilizer for baby-bottle (not the micro-wave type). Perform the process faithfully following the instructions of the sterilizer. To ensure that the disinfection is effective, choose a sterilizer with an operating cycle of at least 6 minutes.

- After disinfection shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).

AIR FILTERING

The unit is equipped with an air filter (A3) which should be replaced when dirty or when its colour changes. To replace the filter take it out by pulling it as shown in figure. The filter is designed to remain fixed in its housing. Use original MANUFACTURER accessories only.

SYMBOLE

Alternating current

CE Medical Marking ref. Dir. 93/42 EEC and subsequent amendments

Important: check the operating instructions

Class II device

"OFF" for part of equipment

"ON" for part of equipment

Serial number of device

In compliance with European standard EN 10993-1 Biological evaluation of medical devices and 93/42/EEC Medical Device Directive. Phthalates free. In conformity with Reg. (EC) no. 1907/2006

Type BF applied part

Hazard: electrocution. Consequence: Death. Do not use device while taking a bath or shower

IP21 Enclosure protection rating: IP21. (Protected against solid bodies over 12 mm. Protected against access with a finger; Protected against vertically falling water drops.)

Hersteller

TECHNICAL DATA

Model: NM-F700

Voltage: 230V~ 50Hz 140VA
Max pressure: 2.6 ± 0.4 bar
Compressor air output: 10 l/min approx
Sound level (at 1 m): Approx. 57 dB (A) with open cover
Operation: approx. 54 dB (A) with closed cover
Operating conditions: continuous
Temperature: min 10°C; max 40°C
RH Air humidity: min 10%; max 95%
Atmospheric pressure: min 69KPa; max 106KPa

Storage conditions: Temperature: min -25°C; max 70°C
RH Air humidity: min 10%; max 95%
Atmospheric pressure: min 69KPa; max 106KPa

Dimensions (W)x(D)x(H) 270x180x110 cm
Weight: 1.7 Kg

APPLIED PARTS


Type BF applied parts are: patient accessories (C2, C3, C4, C5)

RF8 NEBULIZER

Medication minimum capacity ml	Medication maximum capacity ml	(1)Operating pressure bar	(1)Delivery ml/min (approx)		(2)MMAD (µm)		(2)EINATEMBARE FRAKTIONEN < 5 µm (FPF) %	
2	8	1.15	max	max	max	max	max	max
			0.44	0.28	4.12	3.69	59.8	64.2

(1) Values measured according to Manufacturer's internal procedure I29-P07.5 Delivery rate values can change depending on the patient's respiratory capacity.

(2) In vitro characterization certified by TÜV Rheinland LGA Product GmbH - Germany in compliance with the European Standard EN 13544-1, ANNEX CC for nebulizing system. More details are available on request.

DEVICE DISPOSAL

In conformity with Directive 2012/19/EC, the symbol shown on the device to be disposed of indicates that it is considered as waste and is therefore subject to "sorted waste collection". The user must therefore deliver (or have delivered) the above waste to an appropriate waste collection centre provided by the local authorities, or deliver it to the dealer when purchasing a new appliance of the same type. Presorted waste collection and the subsequent treatment, recovery and disposal operations favor the production of appliances made of recycled materials and limit the negative effects to the environment and public health due to incorrect waste management. The unlawful disposal of the product by the user could result in administrative fines as provided by the laws transposing Directive 2012/19/EC of the European member state or of the country in which the product is disposed of.

ELECTROMAGNETIC COMPATIBILITY

This device has been designed to satisfy requirements currently required for electromagnetic compatibility (EN 60 601-1-2:2007). Electrical medical devices require special care during installation and use with respect to EMC requirements. It is therefore required that they be installed and/or used according to the manufacturer's specification. Potential risk of electromagnetic interference with other devices, in particular with other devices for analysis and treatment. Radio and mobile telecommunications devices or portable RF (mobile phones or wireless connections) may interfere with the operation of electrical medical devices. For further information visit www.new-med.at. Manufacturer reserves the right to make technical and functional changes to the product without notice.

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NEWMED GmbH